C-Peptide ELISA Kit

Code K6220

An enzyme immunoassay for the quantitative measurement of C-peptide in human clinical samples. This kit contains reagents for 96 test wells.

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1 INTENDED USE

Dako C-peptide ELISA Kit is an enzyme-linked immunosorbent assay (ELISA) for the quantitative measurement of C-peptide in human serum and plasma.

2 SUMMARY

C-peptide is formed as a by-product of the processing of insulin in the pancreatic β cell. Proinsulin is cleaved proteolytically to form equimolar amounts of mature insulin and C-peptide which are released into the portal circulation. So called because it connects the A and B chains of insulin in the proinsulin molecule, C-peptide is a single chain of 31 amino acids (mol wt 3020 D) and, unlike insulin, has no known physiological function. Because C-peptide has a longer half-life than insulin (2-5 times longer) higher concentrations of C-peptide persist in the peripheral circulation, and these levels fluctuate less than insulin. For these reasons, plasma C-peptide concentrations may reflect pancreatic insulin secretion more reliably than the level of insulin itself.¹

C-peptide measurement is used for the following clinical applications:-¹

1. To assess the residual β cell function in patients treated with insulin and to distinguish between type 1 diabetes and type 2 diabetes. Of particular interest is its use to indicate the need for progression to insulin therapy in type 2 diabetes.
2. The diagnosis of factitious hypoglycaemia. The surreptitious administration of insulin causes high insulin levels in the absence of elevated C-peptide concentrations.
3. Insulinoma diagnosis, especially in patients treated with insulin. C-peptide measurement is used in insulin suppression tests in euglycaemic patients with suspected insulinoma. Elevated C-peptide levels in this test are indicative of insulinoma.
4. As a marker for residual pancreatic tissue after pancreatectomy. In the case of insulinoma, C-peptide measurement may be used to detect metastasis and the response to therapy. It may also be used to monitor the progress of pancreas or islet cell transplantation.


Until recently routine C-peptide measurement has only been possible by radioimmunoassay, although HPLC methods have been described in the literature.² Dako C-peptide is an enzyme immunoassay for the quantitative determination of C-peptide.

3 PRINCIPLE OF THE TEST

Dako C-peptide is an ELISA based on two monoclonal antibodies. Simultaneous incubation of specimen and enzyme-labelled antibody in a microplate microwell coated with a specific anti-C-peptide antibody forms a complex. A simple washing step removes unbound enzyme-labelled antibody. The bound conjugate is detected by reaction with the substrate 3,3',5,5'-tetramethylbenzidine (TMB). The reaction is stopped by adding acid to give a colorimetric endpoint that is read spectrophotometrically.

The inclusion of calibrators of known C-peptide concentration in the assay allows a calibration curve to be constructed from which the level of C-peptide in samples can be determined. The assay principle is shown below in Figure 1.

Figure 1: Dako C-peptide assay principle

4 DEFINITIONS

Explanation of symbols:

REF	Product code and catalogue number
	Consult the instructions for use
	Contains sufficient for 'N' tests
	Manufacturer
IVD	In vitro diagnostic medical device
	Use by
LOT	Batch code
	Storage temperature limitations
	Add water

5 REAGENTS PROVIDED

96 - Each kit contains sufficient materials for 96 determinations. - The shelf life of the kit is as indicated on the outer box label.

5.1 DAKO C-PEPTIDE TEST CONTENTS

One Instructions for Use Booklet.

MICROTITRATION PLATE One x 96 microwell plate of 12, 8 microwell strips coated with mouse monoclonal anti-C-peptide antibody. A resealable plastic pouch is provided for storage of unused microwells.

One bottle of each of the following:

CALIBRATOR 1 5 mL Calibrator 1: Protein matrix with an anti-microbial agent.

CALIBRATOR 2 Calibrators 2 - 5 (lyophilised): Synthetic C-peptide in a protein

(109467-003)

K6220/EFG/CKJ/2009.06.17 p. 3/43

CALIBRATOR 3
CALIBRATOR 4
CALIBRATOR 5

matrix with an anti-microbial agent. The assay is calibrated against International Reference Reagent IRR 84/510, and values are expressed in both ng/mL and in S.I. units (pmol/L) (calculated from the molecular weight of C-peptide, 3020 D). Precise values in both units are stated on the Value sheet included with the product and must always be used for calibration of the assay. The approximate values are:-

Calibrator	1	2	3	4	5
ng/mL	0	0.3	3	9	15
pmol/L	0	100	1000	3000	5000

CONJUGATE CONCENTRATE

2.5 mL Conjugate Concentrate: Anti-C-peptide monoclonal antibody conjugated to horseradish peroxidase (HRP) in an organic buffered solution containing protein, inorganic salts, a red coloured dye and an anti-microbial agent.

CONJUGATE DILUENT

10 mL Conjugate Diluent: Organic buffered solution containing inorganic salts, protein and an anti-microbial agent.

WASH BUFFER (15x)

28 mL Wash Buffer Concentrate (15x): Organic buffered solution containing detergent, inorganic salts and an anti-microbial agent.

SUBSTRATE TMB

13 mL Substrate: stabilised peroxide and 3,3'-5,5'-tetramethylbenzidine (TMB) substrate in a dilute buffer solution. TMB has been reported to be non-carcinogenic. However, personal protective equipment is recommended to avoid direct exposure.

STOP SOLUTION

13 mL Stop Solution: 0.46 mol/L sulphuric acid.

5.2 RECONSTITUTION, STORAGE AND RE-USE OF KIT COMPONENTS

Dako C-peptide allows for testing of up to 86 batches of specimens during any 8 week period. In order to ensure optimal kit performance, it is important that all reconstituted and unused kit components are stored according to the following instructions.

5.2.1 Monoclonal Antibody Coated Microwells - MICROTITRATION PLATE

Open the plate pouch by cutting along the seal. Break off the required number of microwell strips and relocate them into the frame. **Use only complete strips.** Place unused microwells in the resealable plastic pouch with the dessicant, carefully reseal the pouch and store at 2-8 °C. Microwells may be used for up to 8 weeks after initial opening, provided they are stored in this manner.

5.2.2 Calibrator 1 - CALIBRATOR 1

Ready to use. Store unused Calibrator 1 at 2-8 °C. Calibrator 1 may be used for up to 8 weeks at 2-8 °C after initial opening, provided it is stored in this manner.

5.2.3 Calibrators 2 - 5 - CALIBRATOR 2CALIBRATOR 3CALIBRATOR 4CALIBRATOR 5

Reconstitute the lyophilised calibrators with 1 mL of fresh distilled or deionised water. Discard the stoppers after reconstitution. Allow to stand for 10-15 minutes at room temperature, mix gently before use. Aliquot and store reconstituted calibrators upright at 2-8 °C for up to 5 days or -20 °C for up to 8 weeks .

5.2.4 Conjugate Concentrate and Conjugate Diluent – CONJUGATE CONCENTRATECONJUGATE DILUENT

Dilute the conjugate concentrate (2.5 mL) by adding the conjugate diluent (10 mL). This is sufficient for 96 microwells. Mix gently before use. Store working strength Conjugate at 2-8 °C for up to 8 weeks.

5.2.5 Wash Buffer Concentrate - WASH BUFFER (15x)

Provided 15x concentrate. Prepare working strength wash buffer by adding 1 part of wash buffer to 14 parts of fresh deionised or distilled water. For 96 microwells dilute 28 mL of wash buffer concentrate with 392 mL of deionised or distilled water to give a total volume of 420 mL.

Prepare working strength wash buffer as required on the day of use. Store unused concentrate at 2-8 °C.

Do not store unused working strength wash buffer for subsequent use (see Section 8.2.11).

5.2.6 Substrate Solution - SUBSTRATE TMB

Ready to use. Store unused substrate at 2-8°C, protect from light.

5.2.7 Stop Solution - STOP SOLUTION

Ready to use. Store unused reagent at 2-8 °C.

6 ADDITIONAL REAGENTS

6.1 REAGENTS

Fresh deionised or distilled water for reconstitution of Calibrators 2-5 and preparation of working strength wash buffer.

6.2 ACCESSORIES

The following products are intended for use in conjunction with Dako C-peptide. Contact your local Dako subsidiary or distributor for further information.

The Diabetes Immunoassay Control Set (Code S6247) is intended as a tri-level control for use with diabetes immunoassays.

7 EQUIPMENT

The following equipment is required:

Precision pipettes and disposable tips or adjustable multichannel pipettes tips to deliver 25 µL -250 µL.

V-shaped reagent troughs.

Timer.

Clean absorbent paper (on to which microwells can be tapped dry).

Measuring cylinder 500 mL.

Microtitration plate shaker. For information on suitability of plate shakers contact your local Dako subsidiary or distributor.

Automated plate washer (optional) or suitable equipment for washing 8 microwell strips (Section 10.2.4).

Spectrophotometer or EIA plate reader able to read a 96 microwell plate of 8 microwell strips at an absorbance of 450 nm with a reference at 620-650 nm. (Optional, Section 10.3, Reading the Test Results).

Application notes for use on open automated systems are available for this assay. Contact your local Dako subsidiary or distributor.

8 PRECAUTIONS

 - For in vitro diagnostic use. For professional use only.

8.1 SAFETY PRECAUTIONS

- 8.1.1** The stop solution is sulphuric acid (0.46 mol/L). Avoid eye and skin contact by wearing protective clothing and eye protection.
- 8.1.2** Clinical specimens may contain infectious agents eg Hepatitis B virus and HIV. All specimens, calibrators, controls and materials that come into contact with them, should be treated as potentially infectious and handled accordingly.
- 8.1.3** Ensure that cuts, abrasions and skin lesions are properly protected and take precautions to avoid self inoculation or splashing of mucous membranes.
- 8.1.4** As with any product derived from biological sources, proper handling procedures should be used.
- 8.1.5** Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
- 8.1.6** Unused solution should be disposed of according to local, State and Federal regulations.
- 8.1.7** Safety data sheet available for professional user on request.

8.2 TECHNICAL PRECAUTIONS

- 8.2.1** Components must not be used after the expiry date printed on the labels. Do not mix different lots of reagents.
- 8.2.2** Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in Section 5.2.
- 8.2.3** Avoid contamination of reagents.
- 8.2.4** Every effort should be made to ensure that no cross-contamination occurs between wells. It is essential that enzyme-conjugated monoclonal antibody is not allowed to contaminate reagents and equipment. A separate pipette should be reserved for dispensing the conjugate and care should be taken to avoid touching or splashing the rim of the well. Conjugate allowed to dry on to the rim of the well may adversely affect the performance of the test.
- 8.2.5** All equipment including washers and liquid handling devices must be free of microbial contamination and be correctly calibrated and maintained according to the manufacturer's instructions.
- 8.2.6** Do not store unused working strength wash buffer for subsequent use. When not in use wash buffer reservoirs should be rinsed in deionised or distilled water and left to dry.
- 8.2.7** Microwells cannot be re-used.
- 8.2.8** Do not use substrate showing a blue colour prior to its addition to the microwells.
- 8.2.9** Protect substrate from light.
- 8.2.10** It is recommended that samples are run in duplicate until the laboratory is totally familiar with the assay.
- 8.2.11** Test results should be interpreted in conjunction with information available from epidemiological studies, clinical assessment of the patient and other diagnostic procedures.

9 COLLECTION AND PREPARATION OF SPECIMENS

Dako C-peptide has been designed for use with human serum and plasma.

9.1 SERUM

Collect blood by venepuncture into tubes without additives, allow to clot and separate the serum from cells by centrifugation. Store samples for up to 3 days at 2-8 °C or for up to 28 days at – 20°C. Aliquot, if necessary, to avoid repeated freezing and thawing of samples, and do not thaw frozen samples by heating above ambient temperature.

9.2 PLASMA

Collect blood by venepuncture into tubes containing heparin or EDTA as anticoagulant. Centrifuge, separate the plasma fraction and store for up to 3 days at 2-8 °C or for up to 28 days at –20 °C. Aliquot, if necessary, to avoid repeated freezing and thawing of samples, and do not thaw frozen samples by heating above ambient temperature.

10 TEST PROCEDURE

PLEASE REFER TO SECTION 8.2 TECHNICAL PRECAUTIONS BEFORE PERFORMING TEST PROCEDURE.

10.1 GENERAL NOTES

10.1.1 Ensure that, within any run, all the reaction times are the same by adding reagents in the same timed sequence.

10.1.2 Pipette into the reagent trough only as much liquid as necessary for any assay.

Do not return excess reagents to bottles after use.

10.2 ASSAY PROCEDURE

NOTE: The assay procedure requires the use of a microtitration plate shaker. For information on suitability of shakers contact your local Dako subsidiary or distributor.

10.2.1 Calibrator, Specimen and Control Addition

Locate the required number of microtitration strips into the microwell holder, according to the number of determinations required. Pipette 25 µL of calibrator, specimen or control into the microwells. Duplicate microwells should be used for the calibrators. **The position of each patient sample on the microtitration plate must be recorded.**

10.2.2 Conjugate Addition

After addition of all calibrators, specimens or controls, add 100 µL of conjugate into each microwell. Avoid cross contamination between microwells.

10.2.3 First Incubation

Place the microtitration plate on the plate shaker and shake during the incubation at 20-30 °C for 60 minutes.

10.2.4 Washing the Microwells

The microwells should be washed using working strength wash buffer (Section 5.2.5).

The washing technique is critical to the test performance and should be carried out so as to ensure complete filling (with a minimum of **250 µL** of working strength wash buffer) and emptying of the microwells.

Three wash cycles are essential, by either automated or manual washing techniques.

Manual washing

If washing microwells manually, aspirate or shake out the contents of the microwells and using freshly prepared wash buffer, ensure complete filling and emptying of microwells. Between each wash step remove all remaining wash buffer by tapping the inverted microwells on to clean absorbent paper. Repeat the process twice, rotating the plate after each wash to start in the opposite direction (a total of 3 washes). Manual washing efficiency is further ensured if the wash buffer is delivered at an angle so as to produce a vortex in the microwells. After the final wash, the plate should be inverted and tapped on absorbent paper to remove the last traces of wash buffer.

Automated washing

Automated washers should be programmed to complete 3 wash cycles. Soak steps are not required. Washers must be correctly calibrated to ensure complete filling and emptying of microwells between each wash. After the final wash, the plate should be inverted and tapped on absorbent paper to remove the last traces of wash buffer.

10.2.5 Substrate Addition and Second Incubation

Add 100 µL substrate solution to each microwell. Place the microtitration plate on the plate shaker and shake for 10 minutes at 20-30 °C.

Note: If using a plate reader with a maximum absorbance of 2.2 AU, reduce the substrate incubation time to 7 minutes.

10.2.6 Stop the Reaction

Add 100 µL stop solution to each microwell, allowing the same time interval between each strip as in Section 10.2.5. Ensure thorough mixing in the microwells.


10.3 READING THE TEST RESULTS

10.3.1 PHOTOMETRIC READING

The plate must be read photometrically within 30 minutes after addition of the stop solution. **Protect the plate from light until ready to read.** Mix the contents of the microwells and read the absorbance of each microwell using a suitable spectrophotometer or EIA plate reader set at 450 nm. Ensure that the bottoms of the microwells are clean before reading and check that no foreign matter is present in the microwells. The reader should be blanked on air (ie with no plate in the carriage) before the plate is scanned.

Alternatively if the spectrophotometer or EIA plate reader allows for the use of a reference wavelength (at 620 to 650 nm), dual wavelength reading should be performed which will eliminate any potential interference caused by aberrations, such as dirt or marks, on the optical surface of the microwells.

10.4 SUMMARY OF DAKO C-PEPTIDE ASSAY PROCEDURE

11 CALCULATION AND RESULTS

Patient sample results are estimated from a calibration curve constructed either manually or automatically by computer.

11.1 AUTOMATED CALCULATION

Refer to the plate reader manual. Plot the absorbance obtained against the concentration of C-peptide in the calibrators. A four parameter curve fit is recommended to construct the calibration curve from which patient sample results are estimated.


Owing to the nature of the four parameter fit, it may be necessary to enter the value of the calibrator 1 as 0.001. This will ensure that the curve passes through all data points. For further details please contact your local Dako subsidiary or distributor.

11.2 MANUAL CALCULATION

Results can be calculated manually using graph paper. Plot C-peptide concentration on a linear scale on the x axis and the absorbance on a linear scale on the y axis. Construct the calibration curve and use this to determine the concentration of C-peptide in patient samples.

11.3 TYPICAL CALIBRATION CURVE

Figure 1: A typical C-peptide calibration curve (ng/mL)

Note: This graph is for the purpose of illustration only and should not be used to determine results.

11.4 DILUTION OF SAMPLES OF HIGH CONCENTRATION

Samples with concentrations of C-peptide expected to be in excess of the top calibrator should be diluted 1 in 4 (v/v) with Calibrator 1 before performing the assay.

12 PERFORMANCE LIMITATIONS

12.1 INTERFERENCES

Clinical evaluations indicate that samples containing high concentrations of lipid or bilirubin do not interfere in the C-peptide assay.

No interference from rheumatoid factor or heterophilic antibodies was seen.

12.2 HAEMOLYSIS

Grossly haemolysed samples should not be used as lower assay values may be obtained.

12.3 "HIGH DOSE HOOK EFFECT"

In two site immunoassay systems the response curve may invert in the presence of extremely high analyte concentrations. In the C-peptide assay this effect occurs at greater than 120 ng/mL (40,000 pmol/L) C-peptide.

12.4 SERUM AND PLASMA

Both serum and plasma can be used. During clinical evaluation, no clinically significant difference was found in measured C-peptide between matched serum and plasma samples. Both heparin and EDTA anticoagulants may be used.

12.5 PROINSULIN AND PROINSULIN COMPLEXES

Even though circulating C-peptide is in 50 times molar excess compared to proinsulin, samples with elevated proinsulin levels will produce increased values in C-peptide (Section 14.2.3).

The presence of circulating proinsulin/anti-insulin antibody complexes in some patients may also interfere with the assay. In such cases it is suggested that the C-peptide levels are measured after precipitation with polyethylene glycol (PEG)³.

12.6 RENAL INSUFFICIENCY

Because the kidney is the major site for C-peptide metabolism, patients with severe renal insufficiency may have abnormally high circulating C-peptide levels.

12.7 FORENSIC APPLICATIONS

The C-peptide assay is not intended for forensic use and has not been validated with post mortem samples.

13 EXPECTED VALUES

It is recommended that each laboratory establishes its own reference intervals.

Reference intervals can differ considerably between methods. Do not use reference intervals established for other methods.

14 PERFORMANCE CHARACTERISTICS

The data presented in this section was obtained during independent clinical evaluations and/or at Dako.

14.1 PRECISION

14.1.1 Within (intra) assay precision

Three pools of patient samples covering the range of the assay were tested 24 times on a single plate.

Pool	n	Mean	SD (ng/mL)	Mean	SD (pmol/L)	CV (%)
1	24	1.55	0.04	513	13.3	2.6
2	24	3.6	0.12	1190	39.3	3.3
3	24	6.2	0.24	2050	80.0	3.9

14.1.2 Between (inter) assay precision

Three pools of patient samples covering the range of assay were tested in 10 assays.

Pool	n	Mean	SD (ng/mL)	Mean	SD (pmol/L)	CV (%)
1	10	1.77	0.08	586	26.4	4.5
2	10	3.4	0.05	1130	17.0	1.5
3	10	5.1	0.17	1690	55.8	3.3

14.1.3 Detection Limit (Sensitivity)

The limit of detection was assessed in 2 assays by testing 20 replicates of the zero calibrator. The average value obtained corresponds to a signal 2.5 standard deviations from the mean signal of the calibrator 1.

Detection limit: 0.09 ng/mL (30 pmol/L)

14.2 ACCURACY

14.2.1 Linearity

Three patient samples (heparin plasma) with high analyte concentrations were diluted with calibrator 1.

Sample	Dilution Factor	Result (ng/mL)	Observed/Expected (%)
1	Neat	7.6	
	1:2	3.8	100
	1:4	2.0	105
2	Neat	8.0	
	1:2	4.1	103
	1:4	2.1	105
3	Neat	8.1	
	1:2	4.1	101
	1:4	2.1	104

14.2.2 Recovery

Three plasma sample pools with low analyte concentrations were spiked with different concentrations of C-peptide and the recovery calculated as a percentage.

Added C-peptide (ng/mL)	Result (ng/mL)	Recovery (%)
Sample pool 1	2.4	
Sample pool 1 + 1.9	4.3	100
Sample pool 1 + 4.3	6.7	100
Sample pool 1 + 7.8	9.7	94
Sample pool 2	1.9	
Sample pool 2 + 1.9	3.6	90
Sample pool 2 + 4.3	5.9	93
Sample pool 2 + 7.8	9.5	97
Sample pool 3	1.5	
Sample pool 3 + 1.9	3.5	105
Sample pool 3 + 4.3	6.0	105
Sample pool 3 + 7.8	9.3	100

14.2.3 Cross reactivity

The specificity of the assay was assessed by observing cross reaction with structurally related peptides and hormones.

Cross reactivity (%) was calculated from the mass of C-peptide and the mass of cross reactant giving the same assay signal.

Analysis	Cross reactivity (%)
Intact human proinsulin (biosynthetic)	84
32-33 split proinsulin	83
Human insulin	0*

* Insulin at a concentration of 240,000 pmol/L was below the detection limit of the C-peptide assay.

Fasting concentrations of intact and split proinsulins are typically only 1-2% of C-peptide concentrations. Cross-reactivity with these molecules, therefore, is not clinically significant.

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1. **Haibach H, Dix J & Shah J (1987).**
Homicide by insulin administration.
Journal of Forensic Sciences, **32**: 1, 208-216.
2. **Heding LG (1975).**
Radioimmunological determination of human C-peptide in serum.
Diabetologia, **11**: 541-548.
3. **Hanning I, Home PD & Alberti KGMM (1985).**
Measurement of free insulin concentrations: the influence of the timing of extraction of insulin antibodies.
Diabetologia, **28**: 831-835.

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TECHNICAL ADVICE AND CUSTOMER SERVICE

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For all enquiries please contact your local Dako subsidiary or distributor.

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Dako Denmark A/S